

Test Menu

TOPIC	DESCRIPTION
Test Name	Ehrlichia, IgG (IFA)
Other Name (s)	Human monocytic ehrlichiosis,
Analyte(s)	Ehrlichia
Test Code	1710
Lab location	Jacksonville location
Department	Virology
Prior Authorization	Requires prior approval from Regional Epidemiology and notification to the testing lab. Contact local County Health Department to start the process for approval
Required Forms	-Test Requisition Form, DH1847. Medical History needed (i.e., onset date, collection date, travel history, symptoms, and tick bite history).
Specimen Sources	Paired Sera* collected at least 3 weeks apart or single convalescent serum plus appropriate exposure and symptom history.
Supplemental Information- Special Specimen Preparation	Paired Sera*: 1. Collect 3-5 mL of blood into a serum-separator tube (tiger/red - topped tube), 2. Collect 3-5 mL of blood in EDTA (lavender-topped tube) this is for whole blood (with paired serum). 3. Collect 1 mL minimum of clean voided urine in labeled sterile container with an external cap and an internal O-ring.
Minimum Volume	3-5 mL of blood
Storage Conditions	Refrigerate specimens at 2-8°C after collection up to 3 weeks; after 3 weeks freeze at - ≤20°C.
Collection Media	serum-separator tube or vacutainer (tiger/red - topped tube)
Specimen Labeling	-Specimen must be labeled with at least two unique patient identifiers, Ex: Name and DOB. -The collection date and time if submitting multiple specimens. -Information on the specimen must match the requisition.
Packaging and Shipping Instructions and Handling	Specimens must be shipped between (2-8°C) or frozen (≤-20°C) on dry ice. Separate multiple specimens into different bags (preferred).
Test Methodology	Indirect Fluorescent Antibody (IFA)
Turnaround Time	5 - 10 days
Result Indicator	Titer and interpretation or no antibody detected.
Unsatisfactory Specimen	Hemolysis, Unlabeled or mislabeled specimens, insufficient quantity for testing, incorrect collection tube/transport media, grossly contaminated specimen, disparity between ID on sample and paperwork, improper collection, storage or transport of specimen, no test requested, test requested is not performed. If required, the absence of patient history. If required, the lack of patient history that is compatible with test requested. Test order cancelled by provider, broken, or leaked in transit, etc.
Interferences and Limitations	lipemia, hemolysis, microbial growth, Some cross-reactivity with anaplasma and other ehrlichial agents, antibiotic treatment.
Additional Information & Notes	Date of onset, tick exposure, clinical symptoms, and recent travel history is required. The first sample should be taken within the first 2 weeks of illness and the second should be taken 2 to 4 weeks later. Requires prior approval from CHD and notification to the testing lab.
Reference Range	N/A
Reference Lab	CDC if needed
Reflex testing	None