

**NOTES:**

- 1. Screw-on caps are wrapped with Parafilm or adhesive tape.**
- 2. If a computer generated patient label is affixed to the primary receptacle, it is placed directly over the manufacturer's label to allow the sample to be viewed through the receptacle.**
- 3. When shipped by air, the primary or secondary receptacles are able to withstand, without leakage, an internal pressure producing a pressure differential of not less than 95 kPa in the range of – 40°C to + 55°C (– 40°F to + 130°F).**
- 4. All containers provided by the Bureau of Laboratories have been certified and can be used for ground OR air transport**
- 5. Enough absorbent material is used to absorb the entire contents of all primary receptacles.**