



ELISA Results		
WN	SLE	Conclude
Pos	Neg	WN infection
Pos	Equiv	WN infection
Neg	Pos	SLE infection
Equiv	Pos	SLE infection
Pos	Pos	WN infection with cross reactive antibody
Pos	Pos	SLE infection with cross reactive antibody
Pos	Pos	Flavivirus infection-cannot determine if WN or SLE is the cause; additional testing required.
Pos	Equiv	Flavivirus infection-cannot determine if WN or SLE is the cause; additional testing required.
Equiv	Pos	Flavivirus infection-cannot determine if WN or SLE is the cause; additional testing required.
Equiv	Equiv	Possible flavivirus infection-cannot determine if WN or SLE is the cause; additional testing required.
Neg	Neg	Antibody not detected; additional testing required

Red - Require PRNT testing

*Some samples reported as genus positive in PRNT will be retested. Final results may not be reported until wk 7.